

# STANDARD OPERATING PROCEDURE

**Title: Liquid Chromatography, Dionex UltiMate 3000 RSLCnano LC System**

**Version #: 1**

**Author: Hui Zhang Lab**

**Date: 05/08/2014**

## Purpose

The purpose of this document is to describe the liquid chromatography (LC) method for quantitative mass spectrometry-based analyses.

## Scope

This procedure describes the setup of the LC and the method parameters. It is specific to the operation of the Dionex UltiMate 3000 RSLCnano LC system.

## Responsibilities

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

## Equipment

- HPLC: Dionex UltiMate 3000 RSLCnano LC system (Thermo Fisher Scientific; cat. # NCS-3500RS)

## Materials

- Injection loop: 20  $\mu$ L Dionex nanoViper sample loop (Thermo Fisher Scientific; cat. # 6826.2420)
- Column compartment: UltiMate 3000 Binary Rapid Separation Nano Flow Pump with Ternary Loading Pump and Column Compartment (Thermo Fisher Scientific; cat. # NCS-3500RS)
- Autosampler: UltiMate 3000 Thermostatted Rapid Separation Pulled Loop Wellplate Sampler (Thermo Fisher Scientific; cat. # WPS-3000)
- Solvent Degasser: UltiMate 3000 Integrated Solvent and Degasser Rack, 4 Channels (Thermo Fisher Scientific; cat. # SRD-3400)

- Trap Column: 300 µm I.D. x 5 mm packed with Acclaim PepMap 100 5 µm, 100 Å C18 (Thermo Fisher Scientific; cat. # 160454)
- Analytical Column: 75 µm I.D. x 25 cm EASY-Spray column packed with Acclaim PepMap RSLC C18, 2 µm (Thermo Fisher Scientific; cat. # ES802)
- Autosampler vials: 9 mm assembled amber autosampler vial kit (Thermo Fisher Scientific; cat. # C5000-196W)
- Autosampler vial inserts: Polyspring, glass conical insert (Thermo Fisher Scientific; cat. # C4010-630)

## Reagents

- Water: Optima LC/MS-grade (Fisher Scientific; cat. # W6-4)
- Acetonitrile: Optima LC/MS-grade (Fisher Scientific; cat. # A955-4)
- Formic Acid: LC-MS Ultra (Sigma-Aldrich; cat. # 14265)

## Solutions

- Loading pump, mobile phase A: 2% ACN/0.1% formic acid in water
- NanoFlow pump, mobile phase A: 2% ACN/0.1% formic acid in water
- NanoFlow pump, mobile phase B: 0.1% formic acid in 90% ACN

## Procedure

### 1. Autosampler method

Draw Speed =	200 [nl/s]
Draw Delay =	5000 [ms]
Disp Speed =	2000 [nl/s]
Dispense Delay =	2000 [ms]
Waste Speed =	4000 [nl/s]
Wash Speed =	4000 [nl/s]
Loop Wash Factor =	2.000
Sample Height =	2.000 [mm]
Puncture Depth =	8.000 [mm]
Wash Volume =	150.000 [µl]
Rinse Between Reinjections =	No
Low Dispersion Mode =	Off
Inject Mode =	µL pick-up
Inject Volume =	6.000 µL

### 2. Loading pump profile

- a. Flow rate: 5.000 µL/min

- b. Flow gradient: Isocratic; 100% A
  - c. Injection valve switch: 4.000 min
- 3. Gradient method:
  - a. Flow rate: 500 nL/min
  - b. Temperature: 42 °C
  - c. Timetable

<i>Retention (min)</i>	<i>Flow (<math>\mu</math>L/min)</i>	<i>%B</i>
0.00	0.500	4.00
5.00	0.500	5.00
64.00	0.500	40.00
67.00	0.500	90.00
77.00	0.500	90.00
78.00	0.500	4.00
90.00	0.500	4.00

## Referenced Documents

N/A