

STANDARD OPERATING PROCEDURE

Title: Mini-Validation of Repeatability

Version #: 1

Author: Hui Zhang Lab

Date: 05/08/2014

Purpose

The purpose of this document is to describe the characterization of a set of assays based on its repeatability of measurement over 5 days. This is to estimate the performance of the assay measured in a complex sample across several days.

Scope

This procedure addresses the preparation and running of samples for generating the validation samples in accordance with CPTAC Assay Characterization Guidance Experiment #2.

Responsibilities

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

Equipment

- Microcentrifuge
- Vacuum centrifuge

Materials

- Water: Optima LC/MS-grade (Fisher Scientific; cat. # W6-4)
- Acetonitrile: Optima LC/MS-grade (Fisher Scientific; cat. # A955-4)
- Formic Acid: LC-MS Ultra (Sigma-Aldrich; cat. # 14265)
- Methanol: Optima LC/MS-grade (Fisher Scientific; cat. # A456-4)
- Ammonium formate (Sigma-Aldrich; cat. # 70221)

- Ammonium hydroxide (Sigma-Aldrich; cat. # 320145)
- Polysulfoethyl A TopTips 100 – 200 µL (Glygen; cat. # TT2SSA.96)

Reagents

- 10 mM Ammonium formate in 25% ACN, pH 3.0
- 500 mM Ammonium formate in 25% ACN, pH 6.8
- 80:15:5 (vol:vol:vol) Methanol: Water: Ammonium hydroxide
- Stable isotope-labeled standards (SIS)
 - Crude SIS formerly N-glycosylated peptides were synthesized (~60% purity) by Thermo Fisher Scientific (PEPotec SRM peptide library). SIS peptides incorporated a fully atom labeled ^{13}C and ^{15}N isotope at the C-terminal lysine (K) or arginine (R) position of each tryptic peptide, resulting in a mass shift of +8 or +10 Da, respectively. Deamidated asparagine residues corresponding to N-glycosylation sites were synthesized as aspartic acid residues. Peptides were provided in 0.1% TFA/50% ACN and stored at -80 °C until use.
 - A stock SIS mix was prepared at a concentration of 1320 pmol/µL for each peptide, followed by de-salting via strong cation exchange (SCX) as indicated in Procedure 1, below. The peptide recovery following SCX clean-up was assumed to be 50%.
- Matrix
 - A background matrix consisting of trypsin digested serum-derived N-glycopeptides was prepared according to the SOP entitled “Trypsin digestion of serum and enrichment of N-glycopeptides using an automated liquid handler”. This background matrix was used for the preparation of the response curves and for the preparation of the mini-validation of repeatability experiments. Given the limitation of a maximum yield of 20 µg of glycopeptides per preparation of 40 µL serum, N-glycopeptides prepared using multiple DART pipette tips packed with hydrazide resin were combined to achieve the required amount of N-glycopeptides that were needed for the repeatability study.

Procedure

1. SCX de-salting of stock SIS mix
 - a. All centrifugation steps are performed at 2,000 rpm for 1.5 min, unless otherwise specified. The maximum binding capacity of each TopTip is 1 mg.
 - i. Condition Polysulfoethyl A TopTip 2x with 330 µL of Methanol
 - ii. Wash 2x with 330 µL of 10 mM Ammonium Formate in 25% ACN, pH 3.0
 - iii. Wash 2x with 330 µL of 500 mM Ammonium Formate in 25% ACN, pH 6.8
 - iv. Wash 2x with 330 µL of 10 mM Ammonium Formate in 25% ACN, pH 3.0
 - v. Wash 2x with 330 µL of Water
 - vi. Wash 4x with 330 µL of 10 mM Ammonium Formate in 25% ACN, pH 3.0
 - vii. Slowly load acidified sample (pH < 3.0) 2x; centrifuge at 1,100 rpm for 5 min

- viii. Wash 6x with 330 μL of 10 mM Ammonium Formate in 25% ACN, pH 3.0
- ix. Allow TopTip to dry out. Elute sample 2x with 300 μL of 80:15:5 (vol:vol:vol) Methanol: Water: Ammonium Hydroxide; centrifuge at 1,100 rpm for 5 min
- x. Dry eluted sample in a vacuum centrifuge

2. Determination of spike levels and preparation of samples

- a. Peptides were multiplexed according to the LLOQ and linear range determined from the response curves in the CPTAC Assay Characterization Guidance Experiment #1 in order to prepare Validation samples at an appropriate concentration. The three SIS spike levels were as follows:

- i. Low: 0.288 pmol/ μL (1.728 pmol on column)
- ii. Medium: 7.5 pmol/ μL (45 pmol on column)
- iii. High: 12 pmol/ μL (72 pmol on column)

These concentrations were selected to approximate 1.5-3.0x LLOQ (Low), 50-100x LLOQ (Medium), and >100x LLOQ (High).

- b. Dilute matrix to 0.25 $\mu\text{g}/\mu\text{L}$ with 2% ACN/0.2% formic acid.
- c. Prepare dilution mix containing 366.67 μL matrix (0.25 $\mu\text{g}/\mu\text{L}$), 110 μL iRT standard diluted 1:10 and 623.33 μL 0.1 % formic acid (final concentration of matrix in dilution mix is 0.083 $\mu\text{g}/\mu\text{L}$ and final dilution of iRT is 1:100). A total volume of dilution mix should be prepared to permit at least 15 injections of each sample.
- d. Prepare SIS stock by taking 84 μL of SIS (660 pmol/ μL) and drying in a vacuum centrifuge. Resuspend in 369.6 μL dilution mix (prepared in Step 1c, above) for a final SIS stock concentration of 150 pmol/ μL .
- e. Prepare "High" concentration sample by adding 9.6 μL of SIS to 110.4 μL of dilution mix.
- f. Prepare "Medium" concentration sample by adding 6.0 μL of SIS to 114 μL of dilution mix.
- g. Prepare "Low" concentration sample by adding 2.3 μL of 1:10 dilution of SIS to 117.7 μL of dilution mix.
- h. Store samples at -80 $^{\circ}\text{C}$ until LC-PRM MS analysis.

3. Execution of LC-PRM MS analysis

- a. Vortex samples, centrifuge briefly and transfer to autosampler vials. Add sufficient volume to each vial for all of the replicate injections.
- b. Perform LC-PRM MS analysis according to the SOPs entitled "Liquid Chromatography, Dionex UltiMate 3000 RSLCnano LC System" and "Parallel Reaction Monitoring (PRM) Mass Spectrometry, Q-Exactive".

4. Run order

- a. To avoid artificially minimizing variability, the run order of the samples should be randomized. To minimize carryover, 1 wash is inserted after the "Low" and "Medium" samples, and 2 washes are inserted after the "High" samples.

Referenced Documents

- SOP Liquid Chromatography, Dionex UltiMate 3000 RSLCnano LC System
- SOP Parallel Reaction Monitoring (PRM) Mass Spectrometry, Q-Exactive
- SOP Trypsin digestion of serum and enrichment of N-glycopeptides using an automated liquid handler